

Galileo® Med Chair

Side-alternating muscle training while sitting.

In contrast to the standing devices or TiltTable systems, the therapy on the Galileo Med Chair takes place in the sitting position. Important target parameters here are precise pelvic floor activation, back mobilization, trunk stabilization, pain prevention and relaxation. Furthermore, Galileo Med Chair can be used for the treatment of neurological disorders. In particular, patients with difficulty standing up or insecure users are able to activate independently with Galileo Med Chair and also in the familiar environment of the home.



TECHNICAL DATA	
Type number	9N0559 I 2 0 I
Classification	Professional
Medical product	Yes
Certificate	CE0123
Holding possibility	Yes
Ext. Control panel	No
Integr. Control panel	Yes
Remote control	No
Frequency (from/to)	220 Hz
Amplitude (from/to)	0+/-6 mm
Max. Acceleration	9.6 g
Stroke	12 mm
Footplate (I/w)	19 x 12.5 in
Dimensions (I/w/h)	26 x 16 x 20 in
Total weight	62 lbs
Max. Load	440 lbs
Smart Coaching	No
Wobble function	Yes
Option and accessories	

Contact us for your personalized service, training and pricing.

